




DAYS HEAVY DUTY WALKER

\$316.00 - \$320.00 Price range: \$316.00 through \$320.00



SKU: DAYS-HD-BLUE-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Ideal for Bariatric Users looking for a mobility aid for Outdoor Use the **Days Heavy Duty Rollator** is Built for *Optimal Durability, Comfort and Manoeuvrability*.

Constructed out of Heavy Duty Powder Coated Aluminium Frame it can withstand weights of up to 200kg, and is easily foldable for storage and transport.

A wide padded seat with that is 49.5CM wide provides a comfortable seat to sit on and recharge your batteries regardless of whether you are waiting in line for your scripts at the Pharmacy, Shopping or having a walk in the park.

A convenient basket under the seat allows you to take all your essentials such as snacks, water, and medicines with you and still have plenty of room for extra shopping.

FEATURES

- 7.5" wheels
- Adjustable height handles
- Rest seat
- Storage bag
- Loop lockable brakes
- Folds easily for storage & transport

SPECIFICATIONS

- **Max User Weight:** 200kg
- **Seat Height:** 570mm
- **Seat Width:** 495mm
- **Seat Depth:** 355mm
- **Overall Height:** 855 - 975mm
- **Overall Width:** 735mm
- **Overall Depth:** 655mm
- **Backrest:** Curved
- **Wheel:** 7.5"
- **Tyre:** PVC, Black
- **Net Weight:** 8.8kg