




DAYS HIGH BACK CLASSIC DAY CHAIR - UTILITY CHAIR

\$405.00 - \$464.00 Price range: \$405.00 through \$464.00



SKU: ZHI-ZH032

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The High Back Chair features a contoured backrest, padded armrests, vinyl upholstery, adjustable height, and a powder-coated steel frame Day Chair. Available in Teal and Mocha.

KEY FEATURES

- A high back orthopaedic chair that offers adjustability, comfort and support
- Depth Adjustable - Ensures changing user needs are accommodated
- Contoured Backrest - High back design provides increased postural support
- Padded Seat and Armrests - Generous support surfaces for better pressure distribution and sitting comfort
- Padded Headrest Cushion - Provides additional comfort for long periods of sitting
- Sturdy Design - High safe working load of 160kg, yet lightweight for easy manoeuvrability
- Transport Wheels - Allows for the chair to be moved with ease by carers
- Height Adjustable - Easy push button adjustment for comfortable positioning and changing needs
- Non-Slip Tips - Provide stability whilst sitting or transferring
- Powder Coated Steel Frame - Easy to clean, strong and durable

SPECIFICATIONS

Overall Width	640mm
Seat Width	490mm
Seat Depth	480mm

Seat Height	455 - 605mm
Back Height	730mm
Maximum User Weight	160kg

Buy this Days High Back Classic Day Chair, and Utility Chair online today at Medimart as we have NDIS Payment method available.