



## **DAYS HIGH BACK CLASSIC DAY CHAIR - UTILITY CHAIR**



**SKU:** ZHI-ZH032

### **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The High Back Chair features a contoured backrest, padded armrests, vinyl upholstery, adjustable height, and a powder-coated steel frame Day Chair. Available in Teal and Mocha.

### **KEY FEATURES**

- A high back orthopaedic chair that offers adjustability, comfort and support
- Depth Adjustable – Ensures changing user needs are accommodated
- Contoured Backrest – High back design provides increased postural support
- Padded Seat and Armrests – Generous support surfaces for better pressure distribution and sitting comfort
- Padded Headrest Cushion – Provides additional comfort for long periods of sitting
- Sturdy Design – High safe working load of 160kg, yet lightweight for easy manoeuvrability
- Transport Wheels – Allows for the chair to be moved with ease by carers
- Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs
- Non-Slip Tips – Provide stability whilst sitting or transferring
- Powder Coated Steel Frame – Easy to clean, strong and durable

### **SPECIFICATIONS**

<b>Overall Width</b>	640mm
<b>Seat Width</b>	490mm
<b>Seat Depth</b>	480mm

<b>Seat Height</b>	455 - 605mm
<b>Back Height</b>	730mm
<b>Maximum User Weight</b>	160kg

Buy this Days High Back Classic Day Chair, and Utility Chair online today at Medimart as we have NDIS Payment method available.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)