




DAYS LINK FOLDABLE BACKREST SELF PROPELLED WHEELCHAIR



SKU: ANI-DAYWHECWRI-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Days Link Self-Propelled Wheelchair

The Days Link Self-Propelled Wheelchair is designed to provide mobility support for both full-time and part-time wheelchair users. It is suitable for use in personal environments, such as at home and while travelling, as well as in professional settings, including clinical environments. This wheelchair is designed for ease of use, comfort, and safety, making it an excellent choice for users requiring reliable and versatile mobility assistance.

Features

- **Versatile Use:** Suitable for both full-time and part-time wheelchair users, providing reliable mobility support in personal and professional environments.
- **Six Width Options:** Available in six seat widths ranging from 36 cm (14") to 48 cm (19"), ensuring a proper fit for various body sizes.
- **Lightweight Design:** Weighs up to 16.6 kg (36.6 lb), making it easy to manually propel or push with assistance.
- **Puncture-Proof Tyres:** Equipped with puncture-proof tyres for worry-free use on various terrains, enhancing reliability and convenience.
- **Rear Anti-Tippers:** Includes rear anti-tippers for additional safety, particularly on slopes and uneven surfaces.
- **Comfortable Hand rims:** Features comfortable-to-grip hand rims that facilitate easy manual propulsion.
- **Adjustable Backrest:** The height-adjustable backrest can be customised for user comfort and support.
- **Half-Folding Back:** The half-folding back design allows for compact storage and easy transportation.

- **Swing-Back Armrests:** Equipped with swing-back armrests for easier transfers and improved access.
- **Heel Strap:** Includes a heel strap for additional foot support and stability.

Specifications

- **Seat Widths:**
 - 41 cm (16.1"), 46 cm (18.1"), 50.8 cm (20"), 55.8 cm (21.9")
- **Seat Depth:** 43 cm
- **Overall Height:** 98 cm
- **External Depth (With Footrests):** 104 cm
- **External Width (Open):** Seat width + 18 cm
- **Chair Weight:** 17.5 kg
- **Safe Working Load:**
 - 130 kg for 16" and 18"
 - 150 kg for 20" and 22"

The Days Link Self-Propelled Wheelchair offers a combination of comfort, safety, and versatility, making it an ideal choice for users who need reliable mobility support. Its range of features and adjustable components ensure a customised and comfortable experience, whether at home, travelling, or in a clinical setting.

For further information go to [medimart.com.au](https://www.medimart.com.au) or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimart.com.au