




DAYS SEAT WALKER



SKU: DAYS-105-

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Days Seat Walker with Handbrakes & Curved Backrest is a versatile and practical walking aid designed for individuals who require additional support when moving around, whether indoors or outdoors. This walker combines stability with comfort, allowing users to maintain their active lifestyle with confidence.

FEATURES

- **Enhanced Mobility:** The walker is engineered to provide smooth and reliable assistance, featuring large wheels that glide effortlessly over different surfaces.
- **Ergonomic Handbrakes:** Equipped with easy-to-use handbrakes, users can control the walker's speed and ensure immediate stopping when necessary.
- **Comfortable Seating:** The walker includes a built-in padded seat, complemented by a curved backrest that provides users with a comfortable spot to rest whenever needed.
- **Adjustable Height:** The handles of the walker are adjustable, allowing users to set the height according to their personal comfort and ergonomic requirements.
- **Compact Folding Design:** For convenience in storage and transportation, the walker can be folded into a more compact form, easily fitting into car trunks or storage spaces.
- **Durable Frame:** The robust construction ensures longevity and consistent performance, supporting users in their daily activities.
- **Added Safety Features:** Non-slip feet and a sturdy frame design increase the overall safety of the walker, giving users peace of mind.

The Days Seat Walker with Handbrakes and curved Backrest is an excellent solution for those who

value their independence but still need a little help with mobility. Its thoughtful design ensures that users can enjoy their daily routines without the walker becoming a hindrance, proving that practicality can be achieved without sacrificing comfort.

SPECIFICATIONS

- Max User Weight: 125kg
- Seat Height: 565mm
- Seat Width: 350mm
- Seat Depth: 320mm
- Width Between Arms: 460mm
- Backrest Height: 800mm
- Overall Height: 840-940mm
- Overall Width: 620mm
- Overall Depth: 645mm
- Backrest: Straight, folding
- Wheel: 8"
- Tyre: PVC

Buy the Days Rollator walking frame online using the [NDIS payment](#) method as well as AfterPay.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au