




## **DAYS SHOPPER 12 TRANSIT WHEELCHAIR**

**\$411.00**



**SKU:** MFI-70258

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

This transit wheelchair makes it easy for carers and loved ones to provide mobility support, with rear wheel locks, attendant handbrakes and a modern design that can be manoeuvred in tight spaces. Available at a great price, this quality chair includes swing-away leg rests that are removable for your convenience.

The Shopper 12 attendant control wheelchair is ideal for carers and family members looking for a practical option for transporting loved ones and clients to different places. Available at a great price, it is especially convenient for people who are fairly mobile, but struggle with standing or walking for longer periods of time.

Attendant propelled wheelchairs (sometimes called a manual wheelchair or lightweight transit wheelchair) are not like self-propelled wheelchairs, but steered and pushed by a helper. This allows a carer or loved one to easily help someone move from place to place, whether they are out in the community, at the shops or at home.

Lightweight yet durable, this aluminium wheelchair can easily be folded up into a compact form for easy transportation inside your car. It features a folding backrest for easy compact storage, comfortable seating, rear wheel locks and attendant brakes for added safety.

The footplates are height adjustable and flip up when not in use, while the removable footrest can swing away for easy transportation, storage or moving in tight spaces.

## FEATURES

- Folding Backrest
- Rear Wheel Locks
- Padded Armrests
- Safety Lap Belt
- Attendant Handbrakes
- Height Adjustable, Flip-Up Footplates
- Swing-Away, Leg Rests: Removable
- Warranty: One Year Limited

## SPECIFICATIONS

- Weight of Chair- Frame only: 11kg
- Weight of Chair- with Leg Rests: 12.5kg
- Overall Seat Width: 482.5mm
- Seat Width (between arms): 432mm
- Seat Depth: 407mm
- Backrest Height: 457mm
- Leg Rest Height Range\*\*: 356-483mm
- Overall Width (footprint): 559mm
- Overall Depth (footprint): 635mm
- Armrest Height (seat to top of handle): 254mm
- Wheels (solid castors) – Front: 203mm x 25.4mm
- Wheels – Rear: 305mm x 57mm
- Maximum User Weight: 100kg