


## **DAYS SHOPPER 8 TRANSIT WHEELCHAIR**



**\$401.00**



**SKU:** MFI-69397

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The Shopper 8 Attendant Propelled Wheelchair is ideal if you are a Carer and need to transport a loved one to various destinations as it is lightweight, compact and easy to transport. Constructed out of durable aluminium it is ultra lightweight at only 9kg and has an incredible 100kg weight capacity. If safety is important to you then you will love how this wheelchair also comes with a safety lap belt and rear wheel locks.

### **FEATURES**

- Folding Backrest
- Fixed Arm Rests
- Rear Wheel Locks
- Padded Armrests
- Safety Lap Belt
- Attendant Handbrakes
- Height Adjustable, Flip-Up Footplates
- Swing-Away, Leg rests
- Warranty - 12mths Limited warranty.

### **SPECIFICATIONS**

- Weight of Chair - frame only 9kg
- Weight of Chair -with leg rests 105kg
- Overall Seat Width: 432mm
- Seat Width (between arms): 381mm

- Seat Depth: 407mm
- Backrest Height: 432mm
- Leg Rest Height Range: 356 - 483mm
- Armrest Height (seat to top of handles): 254mm
- Wheels (solid castors) Front: 203mm x 25.4mm
- Wheels- Rear: 203mm x 25.4mm
- Maximum User Weight: 100kg