

## DAYS SHOPPER 8 TRANSIT WHEELCHAIR



**SKU:** MFI-69397

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Shopper 8 Attendant Propelled Wheelchair is ideal if you are a Carer and need to transport a loved one to various destinations as it is lightweight, compact and easy to transport. Constructed out of durable aluminium it is ultra lightweight at only 9kg and has an incredible 100kg weight capacity. If safety is important to you then you will love how this wheelchair also comes with a safety lap belt and rear wheel locks.

### FEATURES

- Folding Backrest
- Fixed Arm Rests
- Rear Wheel Locks
- Padded Armrests
- Safety Lap Belt
- Attendant Handbrakes
- Height Adjustable, Flip-Up Footplates
- Swing-Away, Leg rests
- Warranty - 12mths Limited warranty.

### SPECIFICATIONS

- Weight of Chair - frame only 9kg
- Weight of Chair -with leg rests 105kg
- Overall Seat Width: 432mm
- Seat Width (between arms): 381mm

- Seat Depth: 407mm
- Backrest Height: 432mm
- Leg Rest Height Range: 356 - 483mm
- Armrest Height (seat to top of handles): 254mm
- Wheels (solid castors) Front: 203mm x 25.4mm
- Wheels- Rear: 203mm x 25.4mm
- Maximum User Weight: 100kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)