




## **DAYS SHOWER COMMODE ATLANTIC WAVE - TOILET CHAIR**

**\$659.00**



**SKU:** NEW-091563451

### **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

The Atlantic Wave Commode & Shower Chair, 18" Closed Seat is manufactured from a high grade silver painted aluminium frame, which offers a rust free product and high strength. It is easily manoeuvrable and has a weight capacity of up to 160kg making it suitable for bariatric patients. Two sizes available: 18" closed seat width option or 22" open seat width option.

### **FEATURES**

- Easy to clean, detachable polyurethane backrest
- For easier cleaning & access armrests are detachable
- Wipe clean, removable over seat
- Easily accessible, clip on and off , white aperture seat
- Fold-up footrest included
- Durable, easy to clean, plastic coated steel frame
- Use with a square bowl (included), disposable bowl or reusable bowl
- Four wheel brake castors for additional safety

### **SPECIFICATIONS**

- Overall Height: 1020mm
- Overall Width: 520mm
- Overall Depth: 780mm

- Seat Height (Padded): 549mm
- Seat Height (Aperture): 527mm
- Seat Width: 440mm
- Seat Depth: 475mm
- Wheels: 130mm
- Maximum User Weight: 160kg

Buy this Atlantic Wave [Commode Toilet Chair](#) & Shower Chair online at Medimart as we have NDIS Payment method available.