



DAYS SHOWER COMMODE ATLANTIC WAVE - TOILET CHAIR



SKU: NEW-091563451

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Atlantic Wave Commode & Shower Chair, 18" Closed Seat is manufactured from a high grade silver painted aluminium frame, which offers a rust free product and high strength. It is easily manoeuvrable and has a weight capacity of up to 160kg making it suitable for bariatric patients. Two sizes available: 18" closed seat width option or 22" open seat width option.

FEATURES

- Easy to clean, detachable polyurethane backrest
- For easier cleaning & access armrests are detachable
- Wipe clean, removable over seat
- Easily accessible, clip on and off, white aperture seat
- Fold-up footrest included
- Durable, easy to clean, plastic coated steel frame
- Use with a square bowl (included), disposable bowl or reusable bowl
- Four wheel brake castors for additional safety

SPECIFICATIONS

Overall Height: 1020mmOverall Width: 520mmOverall Depth: 780mm



• Seat Height (Padded): 549mm • Seat Height (Aperture): 527mm

• Seat Width: 440mm • Seat Depth: 475mm • Wheels: 130mm

• Maximum User Weight: 160kg

Buy this Atlantic Wave Commode Toilet Chair & Shower Chair online at Medimart as we have NDIS Payment method available.

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au