




## DAYS SUPERLITE DELUXE ROLLATOR

**\$402.00**



**SKU:** HCT-9137D

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Days Superlite Deluxe Rollator is sleek, modern, lightweight, easy to manoeuvre, and folds neatly for an easy functioning rollator. The large front wheels provide added stability and is ideal for outdoor walks. This lightweight, sturdy rollator collapses in one swift motion, folds sideways and stands for easy transport and storage. The design features a storage bag, walking stick/umbrella holder, convenient push-button mechanism which allows for quick height adjustment, lockable loop brakes for security and convenience, and a hidden brake cable design provides added safety. Equipped with a built-in rest seat with a canvas backrest. Available in ruby red or blue.

### FEATURES

- Lightweight, sleek, modern design
- Large front wheels provide better stability. Easy to manoeuvre & ideal for outdoor walks
- Lockable loop brakes for security & convenience
- Side folds neatly, and stands for transport & storage
- Quick height adjustment, hidden brake cable, storage bag
- Equipped with a built-in rest seat with a canvas backrest
- Walking stick or umbrella holder
- Complies with ISO 11199-2:2005

### SPECIFICATIONS

<b>Handle Height</b>	85 - 97cm
<b>Overall Depth</b>	70.5cm

<b>Seat Height</b>	54cm
<b>Backrest Height</b>	29.5 - 42cm
<b>Front wheel diameter</b>	25.4cm
<b>Rear wheel diameter</b>	20.3cm
<b>Seat Width</b>	46cm
<b>Seat Depth</b>	25cm
<b>Product Weight</b>	6.91kg
<b>Maximum User Weight</b>	136kg
<b>Warranty</b>	1 Year Product Warranty, 2 Year Main Frame Warranty
<b>Complies with</b>	ISO 11199-2:2005