



DAYS SUPERLITE DELUXE ROLLATOR

\$402.00 Original price was: \$402.00.<u>\$362.00</u>Current price is: \$362.00.



SKU: HCT-9137D

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Days Superlite Deluxe Rollator is sleek, modern, lightweight, easy to manoeuvre, and folds neatly for an easy functioning rollator. The large front wheels provide added stability and is ideal for outdoor walks. This lightweight, sturdy rollator collapses in one swift motion, folds sideways and stands for easy transport and storage. The design features a storage bag, walking stick/umbrella holder, convenient push-button mechanism which allows for quick height adjustment, lockable loop brakes for security and convenience, and a hidden brake cable design provides added safety. Equipped with a built-in rest seat with a canvas backrest. Available in ruby red or blue.

FEATURES

- Lightweight, sleek, modern design
- Large front wheels provide better stability. Easy to manoeuvre & ideal for outdoor walks
- Lockable loop brakes for security & convenience
- Side folds neatly, and stands for transport & storage
- · Quick height adjustment, hidden brake cable, storage bag
- Equipped with a built-in rest seat with a canvas backrest
- · Walking stick or umbrella holder
- Complies with ISO 11199-2:2005

SPECIFICATIONS

Handle Height	85 - 97cm
Overall Depth	70.5cm



Seat Height	54cm
Backrest Height	29.5 - 42cm
Front wheel diameter	25.4cm
Rear wheel diameter	20.3cm
Seat Width	46cm
Seat Depth	25cm
Product Weight	6.91kg
Maximum User Weight	136kg
Warranty	1 Year Product Warranty, 2 Year Main Frame Warranty
Complies with	ISO 11199-2:2005