




DAYS SWIFT LITE ATTENDANT PROPELLED WHEELCHAIR

~~\$420.00~~ Original price was: ~~\$420.00~~. ~~\$378.00~~ Current price is: \$378.00. Sale Available Online Only!



SKU: ANI-ALTS1703-15

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Days Swift Lite Attendant Propelled Wheelchair is ideal for occasional use. The Swift Lite is designed with the carer in mind, the compact lightweight frame makes it easy to transfer in and out of a car. The push button height adjustable armrests make it versatile allowing you to push under a table then quickly adjust for comfort. Attendant propelled in one size.

KEY FEATURES

- Net Weight 9.5kg
- Folding backrest
- Rear wheel locks
- 17" seat width
- Swing-away, removable leg rests
- Height adjustable, flip-up footplates

SPECIFICATIONS

Days Swift Lite Wheelchair, Attendant Propelled

Overall Height:	895mm
Overall Width:	555mm
Overall Depth:	1050mm
Seat Width:	430mm
Seat Depth:	425mm
Seat Height at Front:	455mm

Seat Height at Rear: 430mm

Buy this Swift Lite Transit Days Wheelchair online with [NDIS Payment Method](#) or Afterpay, shop Now!