



DAYS SWIFT SELF PROPELLED WHEELCHAIR



SKU: ANI-DAYWHE338SP

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Elevate your mobility with the Days Swift Self-Propelled Wheelchair, a versatile and dependable choice for those seeking independence and comfort while navigating various environments. Designed for users who value autonomy and ease of movement, this wheelchair combines robust functionality with user-friendly features.

FEATURES

- **Self-Propelled Wheels:** Equipped with large, sturdy rear wheels, the Days Swift allows users to move themselves independently, providing a sense of freedom and agency.
- **Adjustable Footrests:** For comfort and support, the footrests can be adjusted to suit the user's leg length, ensuring an optimal sitting position.
- **Padded Upholstery:** The seat and backrest include padding to enhance comfort, making extended periods of sitting more enjoyable.
- **Folding Frame:** The wheelchair's frame is designed to fold easily, facilitating transportation and storage in tight spaces.
- **Safety Features:** Comes with integrated brakes for both user and attendant, providing safety and peace of mind during transfers and when stationary.
- **Durable Construction:** Made with a high-quality metal frame, this wheelchair is built to endure daily use and last for years to come.

The Days Swift Self-Propelled Wheelchair is a practical solution for those seeking a wheelchair that supports an active lifestyle. With its combination of self-propulsion, comfort features, and portability, this wheelchair aims to assist users in maintaining independence and mobility in their daily activities.

SPECIFICATIONS

Specifications	35 cm (14") (Paediatric)	40.6 cm (16")	45.72 cm (18")
Overall Height:	91 cm	93 cm	93 cm
Overall Width:	56 cm	62 cm	66 cm
Overall Depth:	98 cm	109 cm	109 cm
Seat Width:	35 cm	41 cm	45 cm
Seat Depth:	40.5 cm	42.5 cm	42.5 cm
Seat Height at Front:	48 cm	50 cm	50 cm
Seat Height at Rear:	45.5 cm	47.5 cm	47.5 cm
Push Handle Height:	91 cm	108 cm	108 cm
Backrest Height from Seat:	40 cm	40 cm	40 cm
Armrest Height from Seat:	22 cm	22 cm	22 cm
Frame Weight:	7.8 kg	7.8 kg	7.8 kg
Leg rest Length:	42 - 54 cm	42 - 54 cm	42 - 54 cm
Front Wheels:	15 cm (6"), solid polyurethane	20 cm (8"), solid polyurethane	20 cm (8"), solid polyurethane
Rear Wheels:	56 cm (22"), solid polyurethane, quick release	61 cm (24"), solid polyurethane, quick release	61 cm (24")
Net Weight:	11.7 kg	12.5 kg	12.5 kg
Maximum User Weight:	80 kg	115 kg	115 kg

For further information go to medimart.com.au or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au