



DAYS SWIFT TRANSIT WHEELCHAIR



SKU: ANI-DAYWHE338T18

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Days Swift Transit Wheelchair is a folding wheelchair designed to be pushed by a carer or attendant, making it ideal for car travel or compact storage when not in use. The push handles have profiled handgrips for a secure and comfortable hold, while its compact dimensions allow for easy manoeuvrability indoors. The ultralight construction and simple folding mechanism make it a versatile chair. The wheelchair has a sturdy aluminium silver frame with 30 centimetre (12 inch) solid rubber rear tyres that are quick release and durable for outdoor use, with no risk of punctures. The folding system is user-friendly, including a fold-down backrest, allowing for easy transport in the car's boot. The chair weighs less than 8 kilograms and has swing-away and removable leg rests with fixed arms. An extended brake handle ensures easy brake application by both the attendant and the user.

FEATURES

- Folding backrest
- Removable leg rests
- Compact for tight spaces
- Swing-away leg rests
- Extended handle brake
- Suitable for indoor and outdoor use

SPECIFICATIONS

- Weight: 7.8kg (frame only) / 12.9kg (total)
- Seat Width: 41cm
- Seat Depth: 42cm

- Backrest Height: 41cm
- Seat Height to Front: 50cm
- Seat Height to Rear: 47.5cm
- Armrest Height: 23cm
- Overall Width: 56cm
- Upholstery: Nylon (flame retardant)
- Weight Capacity: 115kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au