



## **DAYS SWIFT TRANSIT WHEELCHAIR**

**\$469.00**



**SKU: ANI-DAYWHE338T18**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The Days Swift Transit Wheelchair is a folding wheelchair designed to be pushed by a carer or attendant, making it ideal for car travel or compact storage when not in use. The push handles have profiled handgrips for a secure and comfortable hold, while its compact dimensions allow for easy manoeuvrability indoors. The ultralight construction and simple folding mechanism make it a versatile chair. The wheelchair has a sturdy aluminium silver frame with 30 centimetre (12 inch) solid rubber rear tyres that are quick release and durable for outdoor use, with no risk of punctures. The folding system is user-friendly, including a fold-down backrest, allowing for easy transport in the car's boot. The chair weighs less than 8 kilograms and has swing-away and removable leg rests with fixed arms. An extended brake handle ensures easy brake application by both the attendant and the user.

### **FEATURES**

- Folding backrest
- Removable leg rests
- Compact for tight spaces
- Swing-away leg rests
- Extended handle brake
- Suitable for indoor and outdoor use

### **SPECIFICATIONS**

- Weight: 7.8kg (frame only) / 12.9kg (total)
- Seat Width: 41cm
- Seat Depth: 42cm

- Backrest Height: 41cm
- Seat Height to Front: 50cm
- Seat Height to Rear: 47.5cm
- Armrest Height: 23cm
- Overall Width: 56cm
- Upholstery: Nylon (flame retardant)
- Weight Capacity: 115kg