



DAYS T-HANDLE WALKING STICK

\$33.00



SKU: DAYS-20052BR

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

A lightweight, aluminium walking stick featuring a soft ergonomic handle designed for maximum comfort.

The height of the stick can be adjusted to match the height of the user – an important feature when selecting the stick that is right for you. An easy, single button height adjustment mechanism allows you to adjust the telescopic leg. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.

FEATURES

• Supplied with anti-slip rubber tip.

• Height adjustment range: 760mm to 990mm.

• Maximum User Weight: 113kg

Bronze colour