



### **DAYS WHIRL TRANSIT WHEELCHAIR 18"**



**SKU:** KUN-WHIRL45TR

## **PRODUCT DESCRIPTION**

The DAYS Whirl Wheelchair Healthcare Whirl is a versatile wheelchair with a range of features to accommodate different user needs. This sleek and stylish wheelchair is practical and will look good when you are out and about.

### **Key Features**

- Quick-release wheels (for transport)
- Removable armrests
- Swing away and removable leg rests
- Solid tyres (no maintenance required)

### **Specifications**

- Weight: 16.8kg (complete)
- User Weight: 120kg
- Seat Width: 450mm
- Seat Depth: 430mm
- Backrest Height: 440mm
- Seat Height - Front: 500mm
- Seat Height - Rear: 475mm
- Front Castor Wheel: 8" x 1" PU
- Rear Wheel: 24" x 1 3/8" PU
- Overall Width: 640mm
- Upholstery Nylon (Flame Retardant)
- Colour: Silver

Buy the DAYS Whirl Transit Wheelchair using AfterPay or NDIS method Today!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)