



## **DAYS WHIRL TRANSIT WHEELCHAIR 18"**

\$430.00



**SKU:** KUN-WHIRL45TR

## **PRODUCT DESCRIPTION**

The DAYS Whirl Wheelchair Healthcare Whirl is a versatile wheelchair with a range of features to accommodate different user needs. This sleek and stylish wheelchair is practical and will look good when you are out and about.

## **Key Features**

- Quick-release wheels (for transport)
- Removable armrests
- Swing away and removable leg rests
- Solid tyres (no maintenance required)

## **Specifications**

• Weight: 16.8kg (complete)

User Weight: 120kgSeat Width: 450mm

• Seat Depth: 430mm

• Backrest Height: 440mm

Seat Height - Front: 500mmSeat Height - Rear: 475mm

Front Castor Wheel: 8" x 1" PU
Rear Wheel: 24" x 1 3/8" PU

• Overall Width: 640mm

Upholstery Nylon (Flame Retardant)

• Colour: Silver



Buy the DAYS Whirl Transit Wheelchair using AfterPay or NDIS method Today!