



## DECPAC MULTIPURPOSE RAMP

**\$879.00 - \$1,035.00** Price range: \$879.00 through \$1,035.00



**SKU:** DM



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Decpac Ramps are lightweight and easily portable, allowing the user to utilise the ramp whenever necessary and in a wide range of situations. They are manufactured using a combination of uni directional fibreglass and carbon fibre materials, giving the ramp extra strength, without extra weight. The multi-panel design with held together with a patented hinge system which is connected by webbing. It folds quickly into a concertina style. The ramp is easy to use and deploy. For safety the surface has a non-slip, gripping surface as well as pads at the base and top of the ramp. No assembly is required.

### FEATURES

- Portable to allow for vehicle transportation and easy storage
- Lightweight Decpac ramps are very strong yet light enough to make them easily portable
- Compact These ramps fold with a concertina action to allow them to fold up into a much smaller size
- Safe The nonslip surface is suitable for all weather conditions
- Simple to Use This allows the ramp to be deployed in sections for immediate use

The Multipurpose Ramp is 84cm wide and can be purchased in either 87cm or 100cm length.

## COMPARE THE FOUR KINDS OF RAMPS

- **Personal** is ideal for climbing curbs or similar obstacles in and around the home

- [Multipurpose](#) is longer and wider, designed for a wider range of obstacles. This makes it most suitable for travel or commercial environments
- [Senior](#) A large range of options (10 sizes) to cope with obstacles from 22.5cm to 87cm in height
- [Edge Barrier Limiter](#) (EBL) Features edge barriers for additional safety

## SPECIFICATION

<b>Size</b>	87cm - 5.8kg, 100cm - 7.5kg
<b>Ramp width</b>	83.5cm
<b>Max height capacity</b>	14.5 cm and 16.5 cm
<b>Folded width</b>	20 cm
<b>Ramp weight</b>	5.8 kg and 7.5 kg
<b>Safe working load</b>	300kg