



## **DECPAC PERSONAL RAMP**

\$745.00 Original price was: \$745.00.\$671.00 Current price is: \$671.00.



SKU: DP070-F01L05N02-M000

## PRODUCT DESCRIPTION

#### **DECPAC PERSONAL RAMP PURPOSE AND BENEFITS**

The Decpac Ramps are lightweight and easily portable, allowing the user to utilise the ramp whenever necessary and in a wide range of situations. They are manufactured using a combination of uni directional fibreglass and carbon fibre materials, giving the ramp extra strength, without extra weight. The multi-panel design withheld together with a patented hinge system which is connected by webbing. It folds quickly into a concertina style. The ramp is easy to use and deploy. For safety, the surface has a non-slip, gripping surface as well as pads at the base and top of the ramp. No assembly is required.

The 70cm Personal ramp is a portable, lightweight ramp that allows easier access in and around the home or workplace. It is ideal for climbing kerbs or similar obstacles.

#### **FEATURES**

- Portable Allows for vehicle transportation and easy storage
- Lightweight Decapc ramps are very strong yet light enough to make them easily portable
- **Compact** These ramps fold with a concertina action to allow them to fold up into a much smaller size
- Safe The non-slip surface is suitable for all weather conditions
- Simple to Use Allows the ramp to be deployed in sections for immediate use

## COMPARE THE FOUR KINDS OF RAMPS

- Personal is ideal for climbing curbs or similar obstacles in and around the home
- <u>Multipurpose</u> is longer and wider, designed for a wider range of obstacles. This makes is most suitable for travel or commercial environments
- Senior A large range of options (10 sizes) to cope with obstacles from 22.5cm to 87cm in height



• Edge Barrier Limiter (EBL) Features edge barriers for additional safety

# **SPECIFICATIONS**

Length: 70cmWidth: 74.5cm

• Max Height: 17.5cm

• Weight: 4.7kg