




## DECPAC PERSONAL RAMP



**SKU:** DP070-F01L05N02-M000

## PRODUCT DESCRIPTION

### DECPAC PERSONAL RAMP PURPOSE AND BENEFITS

The Decpac Ramps are lightweight and easily portable, allowing the user to utilise the ramp whenever necessary and in a wide range of situations. They are manufactured using a combination of uni directional fibreglass and carbon fibre materials, giving the ramp extra strength, without extra weight. The multi-panel design withheld together with a patented hinge system which is connected by webbing. It folds quickly into a concertina style. The ramp is easy to use and deploy. For safety, the surface has a non-slip, gripping surface as well as pads at the base and top of the ramp. No assembly is required.

The 70cm Personal ramp is a portable, lightweight ramp that allows easier access in and around the home or workplace. It is ideal for climbing kerbs or similar obstacles.

### FEATURES

- **Portable** - Allows for vehicle transportation and easy storage
- **Lightweight** - Decapc ramps are very strong yet light enough to make them easily portable
- **Compact** - These ramps fold with a concertina action to allow them to fold up into a much smaller size
- **Safe** - The non-slip surface is suitable for all weather conditions
- **Simple to Use** - Allows the ramp to be deployed in sections for immediate use

### COMPARE THE FOUR KINDS OF RAMPS

- **Personal** is ideal for climbing curbs or similar obstacles in and around the home
- **Multipurpose** is longer and wider, designed for a wider range of obstacles. This makes is most suitable for travel or commercial environments
- **Senior** A large range of options (10 sizes) to cope with obstacles from 22.5cm to 87cm in height
- **Edge Barrier Limiter** (EBL) Features edge barriers for additional safety

## SPECIFICATIONS

- Length: 70cm
- Width: 74.5cm
- Max Height: 17.5cm
- Weight: 4.7kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)