



DECPAC RAMP WITH EDGE BARRIER



SKU: RAMPEBL

PRODUCT DESCRIPTION

DECPAC EDGE BARRIER RAMP PURPOSE AND BENEFITS

The Decpac Ramp with Edge Barrier is a versatile solution designed to address accessibility challenges for wheelchair and mobility scooter users. With its unique patented hinge system, it combines safety, durability, and convenience, ensuring reliable access to steps, curbs, and other raised surfaces. The edge barrier provides an added layer of protection, ensuring users don't veer off the sides.

Construction & Design:

Crafted from a blend of uni-directional fibreglass and carbon fibre materials, this ramp boasts robust strength without unnecessary weight. Its multi-panel design folds quickly into a concertina style, held together by a patented hinge system connected by webbing.

DECPAC EDGE BARRIER RAMP FEATURES

- **Portable & Compact:** Lightweight construction ensures effortless transport and storage. The concertina folding action reduces its size significantly when stored.
- **Edge Barrier:** Elevated side barriers offer enhanced safety, preventing any accidental roll-offs.
- **Durable Build:** Capable of bearing heavy loads without compromising its shape or integrity.
- **Non-Slip Surface:** A textured grip guarantees secure wheel traction in various weather conditions.
- **Variable Length Options:** Available in widths of 82cm and lengths ranging from 90cm to 350cm.
- **Ready to Use:** No assembly needed – simply unfold and deploy.

Applications:

Perfect for individuals using wheelchairs, mobility scooters, or any mobility device. Whether at home, work, school, healthcare centres or on the go, this ramp is your reliable accessibility companion.

COMPARE THE FOUR KINDS OF RAMPS

- **Personal** is ideal for climbing curbs or similar obstacles in and around the home
- **Multipurpose** is longer and wider, designed for a wider range of obstacles. This makes it most suitable for travel or commercial environments
- **Senior** A large range of options (10 sizes) to cope with obstacles from 22.5cm to 87cm in height
- **Edge Barrier Limiter** (EBL) Features edge barriers for additional safety

SPECIFICATIONS

Length	Width	Weight
90cm	82cm	6kg
120cm	82cm	9.7kg
165cm	82cm	13kg
200cm	82cm	13.7kg
300cm	82cm	27.3kg
350cm	82cm	31.5kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au