



DECPAC SENIOR RAMP

\$852.00 - \$4,557.00 Price range: \$852.00 through \$4,557.00



SKU: RAMPSNR

PRODUCT DESCRIPTION

DECPAC SENIOR RAMP PURPOSE AND BENEFITS

Decpac Senior Ramp is a folding lightweight fibre glass or carbon fibre ramp system. Available in four different styles depending on the needs of the client, these high quality products give users safety and security of functionality.

The Senior is a flexible wheelchair ramp designed to cope with obstacles from 22.5cm – 87cm in height. It is 75cm wide and can purchased with lengths between 90cm - 350cm.

FEATURES

- **Portable** - Allows for vehicle transportation and easy storage
- **Lightweight** - Decapc ramps are very strong yet light enough to make them easily portable
- **Compact** - These ramps fold with a concertina action to allow them to fold up into a much smaller size
- **Safe** - The non-slip surface is suitable for all weather conditions
- **Simple to Use** - Allows the ramp to be deployed in sections for immediate use

COMPARE THE FOUR KINDS OF RAMPS

- **Personal** is ideal for climbing curbs or similar obstacles in and around the home
- **Multipurpose** is longer and wider, designed for a wider range of obstacles. This makes it most suitable for travel or commercial environments
- **Senior** A large range of options (10 sizes) to cope with obstacles from 22.5cm to 87cm in height
- **Edge Barrier Limiter** (EBL) Features edge barriers for additional safety

SPECIFICATIONS

Code	Type	Length	Width	Max Height	Ramp Weight	Max User Weight
DS90	Senior	90cm	75cm	23cm	4.8kg	300kg
DS120	Senior	120cm	75cm	30cm	8.3kg	300kg
DS135	Senior	135cm	75cm	34cm	8.5kg	300kg
DS150	Senior	150cm	75cm	38cm	9.6kg	300kg
DS165	Senior	165cm	75cm	41cm	10.8kg	300kg
DS185	Senior	185cm	75cm	46cm	11.2kg	300kg
DS200	Senior	200cm	75cm	50cm	11.7kg	300kg
DS250	Senior	250cm	75cm	63cm	14.8kg	300kg
DS300	Senior	300cm	75cm	75cm	22.1kg	300kg
DS350	Senior	350cm	75cm	88cm	24.7kg	300kg