


## DECPAC SENIOR RAMP



**SKU:** RAMPSNR

## PRODUCT DESCRIPTION

### DECPAC SENIOR RAMP PURPOSE AND BENEFITS

Decpac Senior Ramp is a folding lightweight fibre glass or carbon fibre ramp system. Available in four different styles depending on the needs of the client, these high quality products give users safety and security of functionality.

The Senior is a flexible wheelchair ramp designed to cope with obstacles from 22.5cm – 87cm in height. It is 75cm wide and can purchased with lengths between 90cm - 350cm.

### FEATURES

- **Portable** - Allows for vehicle transportation and easy storage
- **Lightweight** - Decapc ramps are very strong yet light enough to make them easily portable
- **Compact** - These ramps fold with a concertina action to allow them to fold up into a much smaller size
- **Safe** - The non-slip surface is suitable for all weather conditions
- **Simple to Use** - Allows the ramp to be deployed in sections for immediate use

### COMPARE THE FOUR KINDS OF RAMPS

- **Personal** is ideal for climbing curbs or similar obstacles in and around the home
- **Multipurpose** is longer and wider, designed for a wider range of obstacles. This makes it most suitable for travel or commercial environments
- **Senior** A large range of options (10 sizes) to cope with obstacles from 22.5cm to 87cm in height
- **Edge Barrier Limiter** (EBL) Features edge barriers for additional safety

## SPECIFICATIONS

Code	Type	Length	Width	Max Height	Ramp Weight	Max User Weight
DS90	Senior	90cm	75cm	23cm	4.8kg	300kg
DS120	Senior	120cm	75cm	30cm	8.3kg	300kg
DS135	Senior	135cm	75cm	34cm	8.5kg	300kg
DS150	Senior	150cm	75cm	38cm	9.6kg	300kg
DS165	Senior	165cm	75cm	41cm	10.8kg	300kg
DS185	Senior	185cm	75cm	46cm	11.2kg	300kg
DS200	Senior	200cm	75cm	50cm	11.7kg	300kg
DS250	Senior	250cm	75cm	63cm	14.8kg	300kg
DS300	Senior	300cm	75cm	75cm	22.1kg	300kg
DS350	Senior	350cm	75cm	88cm	24.7kg	300kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)