




DELTA S34 SWIVEL SHOWER STOOL

\$144.00



SKU: 11007

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Swivel Shower Stool makes showering and cleaning easier thanks to its many adjustable and safety features. It comes with padded swivel seat making it possible for the client to turn around without twisting their body.

FEATURES

1. **Height adjustable legs** for tailoring to users requirements.
2. **Padded swivel seat** that makes turning more comfortable and doesn't require twisting of the body.
3. **Suction cup feet** to keep the stool stationary in the shower and reducing the chance of the chair moving when the user initially sits down.
4. **Removable legs** for convenient and compact transportation of the stool.

SPECIFICATIONS

- Product weight: approximately 3kg.
- Width/Depth: 350mm.
- Height: 420-520mm.