


DELUXE BEDSIDE BARIATRIC COMMODE

\$1,055.00 - \$1,173.00 Price range: \$1,055.00 through \$1,173.00



SKU: 1209052-H

PRODUCT DESCRIPTION

Purpose & Benefits of Deluxe Bedside Bariatric Commode

The Deluxe Bedside Bariatric Commode is a highly functional and essential healthcare device designed to assist individuals with mobility challenges, particularly those who need sturdy support due to higher weight considerations. This commode offers a comfortable and secure solution for users, enhancing their independence and dignity in daily personal hygiene activities.

Features

- **Enhanced Comfort:** Comes with padded arms, seat, and backrest, ensuring comfort during use, which is crucial for users spending significant time on the commode.
- **Adjustable Height:** The height of the commode can be adjusted from 475-525mm to fit different user needs, enhancing ergonomic use and accessibility.
- **Durable Construction:** Built to accommodate users up to a certain weight capacity, it is made from high-quality materials that ensure stability and long-lasting use.
- **Easy Accessibility:** Features a hinged seat that allows easy access to the toileting bowl, simplifying the process for users and caregivers.
- **Privacy and Aesthetics:** Equipped with a front and side skirt to discreetly hide the bucket underneath, maintaining the dignity of the user and the aesthetics of the environment.

Specifications

- **Sizes:** GreyStone - 520 mm, Champagne - 600 mm
- **Safe Working Load:** 250 kg
- **Available Colours:** Greystone, Champagne

This Deluxe Bedside Bariatric Commode is particularly valuable in settings such as assisted living facilities, hospitals, and homes where users require reliable support for daily hygiene tasks. Its robust design, coupled with comfort-enhancing features, ensures that it meets the needs of a diverse user base effectively.