



# DEUTSCHER "BIG TED" BARIATRIC FLOORLINE BED - KING SINGLE



SKU: WMR1090BT

#### PRODUCT DESCRIPTION

## **Purpose & Benefits**

The **Deutsche Big Ted Bariatric Bed** is a robust, floor-line electric bed designed for bariatric patients requiring enhanced weight capacity and safety. With a safe working load of 350 kg, it provides superior support while remaining compact and easy to transport or store. Its adjustable height and multiple positioning functions improve patient comfort and assist carers during daily care routines, making it ideal for hospitals, aged care and home care environments.

### **Features**

- Electrically operated 3-function bed high/low, backrest and knee break
- Floor-line design lowers to 135 mm for fall prevention and raises to 735 mm for carer comfort
- Collapsible frame allows easy storage and transport
- Twin actuators on backrest and knee bend sections for added strength
- User-friendly control via hand pendant or integrated foot tower keypad
- Heavy-duty control box with optional battery backup
- Hospital-grade 125 mm dual castors with easy lock/unlock control bars
- Tool-free removal of head and foot boards (Natural Oak Matt finish)



- Air mattress pump bracket fitted to foot end as standard
- Tested to meet ISO 60601-2-52:2010 standards

# **Specifications**

Specification	Details
Model	Deutsche Big Ted Bariatric Bed
Safe Working Load	350 kg (including patient and bedding)
Functions	High/Low, Backrest and Knee Break
Height Range	135 mm - 735 mm
Mattress Size	1070 mm $\times$ 2030 mm (King Single)
Bed Dimensions	2360 mm (L) $\times$ 1070 mm (W)
Bed Weight	115 kg
Castors	125 mm dual hospital-grade wheels with lock control bars
Head & Foot Boards	Full-width removable boards, Natural Oak Matt finish
Compliance	ISO 60601-2-52:2010
Optional Accessories	Side rails (ordered separately), battery backup