



DEUTSCHER "BIG TED" BARIATRIC FLOORLINE BED - KING SINGLE



SKU: WMR1090BT

PRODUCT DESCRIPTION

Purpose & Benefits

The **Deutsche Big Ted Bariatric Bed** is a robust, floor-line electric bed designed for bariatric patients requiring enhanced weight capacity and safety. With a safe working load of 350 kg, it provides superior support while remaining compact and easy to transport or store. Its adjustable height and multiple positioning functions improve patient comfort and assist carers during daily care routines, making it ideal for hospitals, aged care and home care environments.

Features

- Electrically operated 3-function bed – high/low, backrest and knee break
- Floor-line design lowers to 135 mm for fall prevention and raises to 735 mm for carer comfort
- Collapsible frame allows easy storage and transport
- Twin actuators on backrest and knee bend sections for added strength
- User-friendly control via hand pendant or integrated foot tower keypad
- Heavy-duty control box with optional battery backup
- Hospital-grade 125 mm dual castors with easy lock/unlock control bars
- Tool-free removal of head and foot boards (Natural Oak Matt finish)

- Air mattress pump bracket fitted to foot end as standard
- Tested to meet ISO 60601-2-52:2010 standards

Specifications

Specification	Details
Model	Deutsche Big Ted Bariatric Bed
Safe Working Load	350 kg (including patient and bedding)
Functions	High/Low, Backrest and Knee Break
Height Range	135 mm – 735 mm
Mattress Size	1070 mm × 2030 mm (King Single)
Bed Dimensions	2360 mm (L) × 1070 mm (W)
Bed Weight	115 kg
Castors	125 mm dual hospital-grade wheels with lock control bars
Head & Foot Boards	Full-width removable boards, Natural Oak Matt finish
Compliance	ISO 60601-2-52:2010
Optional Accessories	Side rails (ordered separately), battery backup

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au