



DEUTSCHER BIG TED BARIATRIC FLOORLINE BED



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Deutscher Big Ted Bariatric Floorline Bed is a specialised medical bed designed to provide maximum comfort and safety for bariatric patients. It combines robust construction with advanced features to accommodate the specific needs of individuals with higher body weight and limited mobility.

FEATURES

- **Low Floorline Design:** The bed can be lowered close to the floor level, minimising the risk of injury from falls, which is particularly important for patients with mobility issues.
- **High Weight Capacity:** Engineered to support a higher weight capacity, making it suitable for bariatric patients.
- **Adjustable Height and Positioning:** Features electrically adjustable height, backrest, and leg rest to provide optimal comfort and support, as well as to assist with patient care and transfers.
- **Sturdy Construction:** Built with durable materials to ensure stability and longevity, even under frequent and heavy use.
- **Easy-to-Use Controls:** Comes with a user-friendly control panel, allowing both caregivers and patients to adjust the bed settings with ease.
- **Safety Rails:** Includes integrated safety rails that can be raised or lowered as needed, providing additional security for the patient.
- **Lockable Casters:** Equipped with lockable casters to ensure the bed remains stationary during use and to facilitate safe and easy movement when required.

The Deutscher Big Ted Bariatric Floorline Bed represents a significant advancement in medical bed design, offering bariatric patients a safe, comfortable, and dignified healthcare experience. Its combination of bariatric-specific features and user-friendly design makes it a crucial asset in healthcare settings, improving the standard of care for patients with unique needs.

SPECIFICATIONS

- Safe working load: 350kg including client mattress & accessories
- Mattress height range: 135 - 735mm
- Mattress platform: Four panel - 1070 x 2030mm
- Total length of bed: 2360mm
- Width of bed: 1070mm
- Mass of bed: 115kg
- Electrical system: Dewert
- Finish: Powder coated - Pearl White
- Head & foot boards: Full-width boards in Natural Oak Matt are standard

[View the other bed sizes here](#)