



DEUTSCHER CONCORD BED



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This 'state of the art' bed meets the requirements of both client and carer for comfort, safety, and reliability. The automatic two speed lift mechanism allows the client to be raised from floor level to the desired height with ease. High quality and smooth operation is achieved by the use of bearings on all moving parts. The mattress base is available in both standard single and king single sizes. Large 100mm diameter dual castors are standard on all Concord beds. These features make the Concord bed ideal for all healthcare environments. The Concord bed is designed and manufactured in Australia from the highest quality materials. It is tested to meet Australian and international standards (IEC 60601-2-52) and is DHHS compliant.

Deutscher Healthcare is a manufacturing company based in Ballarat, Victoria and has been in continuous operation for over 114 years. In earlier years, the parent company, HF. Deutscher Pty. Ltd., specialised in farm equipment and later moved into producing a range of mowers and outdoor power equipment, some of which is still manufactured today. The company is still owned and operated by the original family.

All products are designed and manufactured to a high quality, using first grade materials to ensure the product you buy meets Australian standard requirements; is safe, reliable and built to give trouble free service for many years to come.

SPECIFICATIONS

MODEL: WMT8080

STD Single: 920mm x2000mm

MODEL: WMT8090

King Single: 1070mm x 2030mm

SAFE WORKING LOAD: 250kg

MAX CLIENT WEIGHT: 210kg

HEIGHT: Maximum raised 810mm / Minimum lowered 135mm

WARRANTY: Frame 10 years / Electrics 3 years

OPTIONS

Battery backup

Underbed Light

Attendant control handset

Different colours & styles of head & foot boards

Attendant control “system” handset

CPR function Handset

Blue tooth Handset

ACCESSORIES

Self-help pole

Bed Sticks

Side rails

Lexia style bed surround package (Fabric covered head & foot board & side rails

Over bed tables

Bed side lockers

Transporter kit