



## DEUTSCHER MENTAL HEALTH BED



**SKU:** N/A

## PRODUCT DESCRIPTION

Available in a single and king single variation.

The Deutscher Mental health bed has been developed to meet the requirements of standard and acute mental healthcare facilities.

- Is raised & lowered when a battery pack is attached.

The bed is strong and robust in design.

Both models are single function bed (Raised & Lowered) to eliminate ligature points.

A generous height range from 240 – 740mm is offered.

The head & foot boards are not removable once the bed is installed.

The bed is available on contract with Hospital Purchasing Victoria.

These beds have been successfully installed in many facilities around Australia.

All electrical wires are encased in a strong metal tube. A stand alone battery pack with hand control is required when raising or lowering the bed.