




DEUTSCHER MENTAL HEALTH BED



SKU: N/A

PRODUCT DESCRIPTION

Available in a single and king single variation.

The Deutscher Mental health bed has been developed to meet the requirements of standard and acute mental healthcare facilities.

- Is raised & lowered when a battery pack is attached.

The bed is strong and robust in design.

Both models are single function bed (Raised & Lowered) to eliminate ligature points.

A generous height range from 240 – 740mm is offered.

The head & foot boards are not removable once the bed is installed.

The bed is available on contract with Hospital Purchasing Victoria.

These beds have been successfully installed in many facilities around Australia.

All electrical wires are encased in a strong metal tube. A stand alone battery pack with hand control is required when raising or lowering the bed.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au