



DEUTSCHER ROSE II BED RAILS - 4 BAR

\$976.00 - \$1,134.00 Price range: \$976.00 through \$1,134.00



SKU: WLR1055S-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The side rails are easily fitted to the beds by two industrial strength clamps. All Medimart side rails are tested to withstand the forces and meet all requirements as of - AS3200.2.38 - 2007.

Raising the latch, at the foot of the bed, allows the rails to fold down. Where an air mattress is to be used, our rails consist of a 4 pole height allowing for maximum safety whilst ensuring the a high level of customisation and comfort. Of course meeting the minimum height requirement.

Side rails are used to minimise the risk of falls from a bed. The side rails can be lowered out of the way when not in use and raised when required.

For Deutscher Beds ONLY.