

DIGITAL BLOOD PRESSURE MONITOR

\$84.00



SKU: UA-704

PRODUCT DESCRIPTION

FEATURES & BENEFITS

Manual inflation Blood Pressure, simple and accurate, low cost, reliable, and our most compact device. The UA-704 is an innovative Palm Top style that makes monitoring your blood pressure easy. The ergonomically designed monitor fits in your palm so you can start your BP measurement easily. Its great portability enables you to carry this unit anywhere to monitor at anytime.

FEATURES

- Innovative Palm Top Style
- Powered by 1 x AA batteries
- Irregular Heartbeat (IHB) Indicator
- SlimFit Comfortable Cuff
- Accuracy Pressure: ± 3 mmHg or 2%, whichever is greater
- Accuracy Pulse: 5%
- Clinically Validated