




## DRIVE AUTOMATIC FOLDING SCOOTER

**\$6,607.00**



**SKU: S3051-H**

## PRODUCT DESCRIPTION

## PURPOSE AND BENEFITS

The latest advanced automatic folding travel mobility scooter from the Solax stable is top of the class. Solax Ezy Auto Fold features remote control folding function, more leg room, adjustable tiller height and angle, comfortable seat and armrests, lithium battery and lightweight!

## FEATURES

- Automatic folding at the press of a button
- Manual folding option
- Remote control
- Lightweight
- Lithium battery
- Easy Operation
- Adjustable tiller height and angle
- Comfortable seat
- Comes with armrests
- Latest Solax technology
- Speed control
- Automatic braking

## SPECIFICATIONS

- **Folding Operation:** Automatic With Remote Control
- **Overall Length:** 930mm

- **Overall Width:** 450mm
- **Tyres:** Solid front 6" Rear 7"
- **Maximum speed:** 6km/hr
- **Speed Control:** Yes
- **Safe Gradient:** 0-12 degrees
- **Weight Capacity:** 125kg
- **Range:** 15km
- **Turning Radius:** 1.4m
- **Weight:** 23.4kg
- **Drive System:** Rear Wheel Drive via Transaxle
- **Seat Height From Ground:** 560mm
- **Battery:** Lightweight Lithium 24v 10ah
- **Motor:** 120w
- **Battery Charger:** 24v 2ah
- **Colours:** Red, Yellow, Blue

Buy the Auto Folding Mobility Scooter Australia online at MediMart using Afterpay or [NDIS method](#) today!