



## **DRIVE POWERSTROLL - DUAL WHEEL REVERSE WHEELCHAIR**

**\$1,920.00 - \$2,155.00** Price range: \$1,920.00 through \$2,155.00



**SKU:** PWCPP009AU-H

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits of Drive Powerstroll - Dual Wheel Reverse Wheelchair**

The **Drive Powerstroll - Dual Wheel Reverse** is an easy-to-use power add-on designed to assist carers in pushing a manual wheelchair with minimal effort. Perfect for use on hills, slopes, or over long distances, the Powerstroll helps reduce fatigue and physical strain, improving safety and comfort for both the carer and the user.

This dual-wheel reverse model features improved traction and power, along with a reverse function for added control in tight spaces. Designed to fit most standard wheelchairs, it quickly attaches and detaches, making it ideal for home, hospital, and community use.

### **Key Features:**

- Motorised drive unit assists carers in pushing manual wheelchairs
- Dual-wheel design for improved traction and performance
- Reverse function for easier manoeuvring in tight spaces
- Fits most standard wheelchairs with quick and simple attachment
- Variable speed control with easy-to-use handset

- Lightweight and portable with convenient battery pack
- Reduces carer fatigue, especially on inclines and longer journeys
- Ideal for home care, hospitals, and aged care environments

## Specifications:

Maximum Gradient:	8°
Motor Size:	24v 150w
Power pack Weight:	19.2kg
Battery Weight:	9kg
Weight without Battery:	10.2kg
Battery Capacity:	12V 12ah
Wheel Size:	200 x 50mm
Speed:	6.5km/h
Range:	16km
Weight Capacity:	135kg