



# DRIVE POWERSTROLL - DUAL WHEEL REVERSE WHEELCHAIR



SKU: PWCPP009AU-H

#### **PRODUCT DESCRIPTION**

# Purpose & Benefits of Drive Powerstroll - Dual Wheel Reverse Wheelchair

The **Drive Powerstroll - Dual Wheel Reverse** is an easy-to-use power add-on designed to assist carers in pushing a manual wheelchair with minimal effort. Perfect for use on hills, slopes, or over long distances, the Powerstroll helps reduce fatigue and physical strain, improving safety and comfort for both the carer and the user.

This dual-wheel reverse model features improved traction and power, along with a reverse function for added control in tight spaces. Designed to fit most standard wheelchairs, it quickly attaches and detaches, making it ideal for home, hospital, and community use.

### **Key Features:**

- Motorised drive unit assists carers in pushing manual wheelchairs
- Dual-wheel design for improved traction and performance
- Reverse function for easier manoeuvring in tight spaces
- Fits most standard wheelchairs with quick and simple attachment
- Variable speed control with easy-to-use handset



- Lightweight and portable with convenient battery pack
- Reduces carer fatigue, especially on inclines and longer journeys
- Ideal for home care, hospitals, and aged care environments

## **Specifications:**

Maximum Gradient:	8°
Motor Size:	24v 150w
Power pack Weight:	19.2kg
Battery Weight:	9kg
Weight without Battery:	10.2kg
Battery Capacity:	12V 12ah
Wheel Size:	200 x 50mm
Speed:	6.5km/h
Range:	16km
Weight Capacity:	135kg

For further information go to medimart.com.au or give us a call: 1300 633 063

## SA

## **Brighton Showroom**

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

#### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

## **Reynella Showroom**

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

#### **VIC**

#### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

#### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au



# NT

## **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 <a href="mailto:support@medimarnt.com.au">support@medimarnt.com.au</a>