



DRIVE SENTRA HEAVY DUTY BARIATRIC SELF-PROPELLED WHEELCHAIR



SKU: MWS442910-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Sentra Bariatric Wheelchair has been designed for the larger person to ensure their mobility experience is as comfortable as possible.

FEATURES

- Comfortable: Comfort is an essential element of its design.
- **Adjustable**: Wheelchair users can easily adjust their seat height to hemi-level with the unique dual-axle position of this chair. The solid castors at the front can also move between 3 different positions for optimal comfort and support.
- Generous Weight Capacity: Up to 200kg (& 318kg for extra wide seat widths)
- **Practical**: Armrests featuring comfortable padding have a comparatively short length that allows access to most tabletops and other surfaces.
- Built-to-Last: Sturdy and stylish wheelchair that has been built-to-last
- **Stylish:** Double-embossed vinyl upholstery and padded desk-style armrests add an attractive look to this chair.
- Low- Maintenance: This wheelchair's design limits tiresome maintenance requirements. A chrome finish added onto the carbon steel frame makes the metal elements of this wheelchair chip-free and easily-maintained. Additionally, the Drive Sentra EC wheelchair is fitted with sealed, precision bearings for a maintenance-free design.
- **Easy Storage & Transportation**: Easy-to-fold to make storage and transportation quick and simple. Additionally, the footrests and armrests can be quickly and easily detached from the rest of the chair for tight storage situations or even sideways transfer of users



SPECIFICATIONS

Code	Display Name			Backrest Height	Seat To Floor Height	I nan	Unit Weight	Colour Material
MWS442910	Heavy Duty Wheelchair - 20in Wide	450 mm	51cm (20")	400 mm	500 mm	125 kg	27.2 kg	Chrome Vinyl
STD22ECDDA-SFAU	Heavy Duty Wheelchair - 22in Wide	450 mm	56cm (22")	400 mm	500 mm	204 kg	27.2 kg	Chrome Vinyl
STD24ECDDA-SFAU	Heavy Duty Wheelchair - 24in Wide	450 mm	61cm (24")	400 mm	500 mm	204 kg	27.2 kg	Chrome Vinyl
STD26DFA-SFAU	Super Heavy Duty Wheelchair - 25.9in Wide	500 mm	66cm (25.9")	450 mm	500 mm	318 kg	36 kg	Chrome Nylon
STD28DFA-SFAU	Super Heavy Duty Wheelchair - 27.9in Wide	500 mm	71cm (27.9")	450 mm	500 mm	318 kg	36 kg	Chrome Nylon
STD30DFA-SFAU	Super Heavy Duty Wheelchair - 29.9in Wide	500 mm	76cm (29.9")	450 mm	500 mm	318 kg	36 kg	Chrome Nylon



For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au