




DRIVE SENTRA HEAVY DUTY BARIATRIC SELF-PROPELLED WHEELCHAIR

\$639.00 - \$3,046.00 Price range: \$639.00 through \$3,046.00



SKU: MWS442910-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Sentra Bariatric Wheelchair has been designed for the larger person to ensure their mobility experience is as comfortable as possible.

FEATURES

- **Comfortable:** Comfort is an essential element of its design.
- **Adjustable:** Wheelchair users can easily adjust their seat height to hemi-level with the unique dual-axle position of this chair. The solid castors at the front can also move between 3 different positions for optimal comfort and support.
- **Generous Weight Capacity:** Up to 200kg (& 318kg for extra wide seat widths)
- **Practical:** Armrests featuring comfortable padding have a comparatively short length that allows access to most tabletops and other surfaces.
- **Built-to-Last:** Sturdy and stylish wheelchair that has been built-to-last
- **Stylish:** Double-embossed vinyl upholstery and padded desk-style armrests add an attractive look to this chair.
- **Low- Maintenance:** This wheelchair's design limits tiresome maintenance requirements. A chrome finish added onto the carbon steel frame makes the metal elements of this wheelchair chip-free and easily-maintained. Additionally, the Drive Sentra EC wheelchair is fitted with sealed, precision bearings for a maintenance-free design.
- **Easy Storage & Transportation:** Easy-to-fold to make storage and transportation quick and simple. Additionally, the footrests and armrests can be quickly and easily detached from the rest of the chair for tight storage situations or even sideways transfer of users

SPECIFICATIONS

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Unit Weight	Colour	Material
MWS442910	Heavy Duty Wheelchair - 20in Wide	450 mm	51cm (20")	400 mm	500 mm	125 kg	27.2 kg	Chrome	Vinyl
STD22ECDDA-SFAU	Heavy Duty Wheelchair - 22in Wide	450 mm	56cm (22")	400 mm	500 mm	204 kg	27.2 kg	Chrome	Vinyl
STD24ECDDA-SFAU	Heavy Duty Wheelchair - 24in Wide	450 mm	61cm (24")	400 mm	500 mm	204 kg	27.2 kg	Chrome	Vinyl
STD26DFA-SFAU	Super Heavy Duty Wheelchair - 25.9in Wide	500 mm	66cm (25.9")	450 mm	500 mm	318 kg	36 kg	Chrome	Nylon
STD28DFA-SFAU	Super Heavy Duty Wheelchair - 27.9in Wide	500 mm	71cm (27.9")	450 mm	500 mm	318 kg	36 kg	Chrome	Nylon
STD30DFA-SFAU	Super Heavy Duty Wheelchair - 29.9in Wide	500 mm	76cm (29.9")	450 mm	500 mm	318 kg	36 kg	Chrome	Nylon

