


DRIVE SENTRA DELUXE HEAVY-DUTY WHEELCHAIR - 26"



SKU: 452144779



PRODUCT DESCRIPTION

Purpose & Benefits

The **Drive Sentra Deluxe Heavy-Duty Bariatric Wheelchair - 26"** is engineered to provide strength, comfort, and independence for larger users who require an extra-wide seat. Its durable carbon steel frame with chrome finish ensures long-lasting performance, supporting users up to 318 kg while maintaining a smooth and stable ride.

Features

- Carbon steel frame with durable chrome coating for a chip-resistant, low-maintenance finish
- Supports a maximum user weight of 318 kg
- Self-propelling design with large rear wheels for comfortable indoor and outdoor use
- Stylish black moulded composite mag wheels with solid tyres for durability
- Plastic hand rims for easy grip and manoeuvrability
- Double embossed vinyl upholstery for comfort and easy cleaning
- Simple folding mechanism for convenient storage and transport

- Height adjustable, swing-away, flip-up steel footrests with ankle straps
- Removable padded armrests for comfort and side transfer ease
- Maintenance-free sealed precision bearings for smooth operation
- Optional elevating footrests available
- Not suitable for use in a motor vehicle

Specifications

	Specification	Value
Seat Width:		66 cm (26")
Seat Depth:		50 cm
Seat Height (min-max):		43 – 48 cm
Overall Width (open):		89 cm
Width (folded):		37 cm
Unit Length:		119 cm
Unit Height:		97 cm
Unit Weight:		40 kg
Safe Working Load (SWL):		318 kg
Warranty:		12 months

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au