



DRIVE SENTRA DELUXE HEAVY-DUTY WHEELCHAIR - 26"

\$2,398.00



SKU: 452144779

PRODUCT DESCRIPTION

Purpose & Benefits

The **Drive Sentra Deluxe Heavy-Duty Bariatric Wheelchair - 26"** is engineered to provide strength, comfort, and independence for larger users who require an extra-wide seat. Its durable carbon steel frame with chrome finish ensures long-lasting performance, supporting users up to 318 kg while maintaining a smooth and stable ride.

Features

- Carbon steel frame with durable chrome coating for a chip-resistant, low-maintenance finish
- Supports a maximum user weight of 318 kg
- Self-propelling design with large rear wheels for comfortable indoor and outdoor use
- Stylish black moulded composite mag wheels with solid tyres for durability
- Plastic hand rims for easy grip and manoeuvrability
- Double embossed vinyl upholstery for comfort and easy cleaning
- Simple folding mechanism for convenient storage and transport



- Height adjustable, swing-away, flip-up steel footrests with ankle straps
- Removable padded armrests for comfort and side transfer ease
- Maintenance-free sealed precision bearings for smooth operation
- Optional elevating footrests available
- Not suitable for use in a motor vehicle

Specifications

Specification	Value
Seat Width:	66 cm (26")
Seat Depth:	50 cm
Seat Height (min-max):	43 – 48 cm
Overall Width (open):	89 cm
Width (folded):	37 cm
Unit Length:	119 cm
Unit Height:	97 cm
Unit Weight:	40 kg
Safe Working Load (SWL):	318 kg
Warranty:	12 months