



MANGAR EAGLE WITH AIRFLO 24

\$6,829.00



SKU: LSS403800

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Mangar Eagle with Airflo 24 is a portable, inflatable lifting cushion that helps people who have fallen get back up to a sitting position. It is ideal for use in homes and other healthcare settings and can be used by a single person to lift a person weighing up to 190kg.

The Eagle with Airflo 24 is a safe and effective way to lift a fallen person. It reduces the risk of back injury for the caregiver and helps to maintain the patient's dignity and independence.

FEATURES

- **High Weight Capacity:** The Eagle with Airflo 24 can lift people weighing up to 190kg, making it a versatile solution for a wide range of users.
- **Airflo 24 Compressor:** The Airflo 24 compressor inflates and deflates the cushion quickly and efficiently, making it easy to use in emergencies.
- **Supportive Design:** The Eagle with Airflo 24 has a built-in backrest that provides additional support to the person being lifted, keeping them comfortable and secure.
- **Portability:** The Eagle with Airflo 24 is lightweight and portable, so it can be easily transported and used in different settings.
- **Single-handed Care:** The Eagle with Airflo 24 can be operated by a single caregiver, making it ideal for use in situations where there is limited assistance available.
- **Risk Reduction:** The Eagle with Airflo 24 reduces the risk of back injury for caregivers by providing a mechanized lifting solution.
- **Maintaining Dignity and Independence:** The Eagle with Airflo 24 is a discreet and effective lifting solution that helps people who have fallen maintain their dignity and independence.

The Mangar Eagle with Airflo 24 is a revolutionary lifting solution that caters to the needs of both caregivers and those they assist. It is a safe, effective, and easy-to-use product that helps to reduce the risk of injury, maintain dignity, and promote independence.

SPECIFICATIONS

- **Maximum lift capacity:** 220 kilograms
- **Inflated seat height:** 56 centimeters
- **Inflated width:** 64 centimeters
- **Back/headrest length:** 52 centimeters
- **Weight of product:** 6.3 kilograms
- **Weight of compressor:** 2.3 kilograms
- **Number of lifts per charge:** Up to 6