




MANGAR EAGLE WITH AIRFLO 24



SKU: LSS403800

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Mangar Eagle with Airflo 24 is a portable, inflatable lifting cushion that helps people who have fallen get back up to a sitting position. It is ideal for use in homes and other healthcare settings and can be used by a single person to lift a person weighing up to 190kg.

The Eagle with Airflo 24 is a safe and effective way to lift a fallen person. It reduces the risk of back injury for the caregiver and helps to maintain the patient's dignity and independence.

FEATURES

- **High Weight Capacity:** The Eagle with Airflo 24 can lift people weighing up to 190kg, making it a versatile solution for a wide range of users.
- **Airflo 24 Compressor:** The Airflo 24 compressor inflates and deflates the cushion quickly and efficiently, making it easy to use in emergencies.
- **Supportive Design:** The Eagle with Airflo 24 has a built-in backrest that provides additional support to the person being lifted, keeping them comfortable and secure.
- **Portability:** The Eagle with Airflo 24 is lightweight and portable, so it can be easily transported and used in different settings.
- **Single-handed Care:** The Eagle with Airflo 24 can be operated by a single caregiver, making it ideal for use in situations where there is limited assistance available.
- **Risk Reduction:** The Eagle with Airflo 24 reduces the risk of back injury for caregivers by providing a mechanized lifting solution.
- **Maintaining Dignity and Independence:** The Eagle with Airflo 24 is a discreet and effective lifting solution that helps people who have fallen maintain their dignity and independence.

The Mangar Eagle with Airflo 24 is a revolutionary lifting solution that caters to the needs of both caregivers and those they assist. It is a safe, effective, and easy-to-use product that helps to reduce the risk of injury, maintain dignity, and promote independence.

SPECIFICATIONS

- **Maximum lift capacity:** 220 kilograms
- **Inflated seat height:** 56 centimeters
- **Inflated width:** 64 centimeters
- **Back/headrest length:** 52 centimeters
- **Weight of product:** 6.3 kilograms
- **Weight of compressor:** 2.3 kilograms
- **Number of lifts per charge:** Up to 6

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

