



ECHO DUAL FUNCTION INFRARED THERMOMETER

\$75.00



SKU: ES-636

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Echo Dual Function Infrared Thermometer is an advanced, non-contact thermometer designed for rapid temperature measurement. It's a versatile tool, suitable for both clinical and household settings, providing accurate readings without the need for physical touch.

FEATURES

- **Dual Functionality:** Measures both body and ambient temperatures.
- **Infrared Technology:** Allows for non-contact measurement, reducing the risk of cross-contamination.
- **Fast and Accurate:** Provides temperature readings within seconds.
- **LCD Display:** The digital screen displays readings clearly and can often have a color-coded fever alarm.
- **Memory Storage:** Capable of storing past temperature readings, which can be beneficial for tracking changes over time.
- **Auto Shut-off:** Conserves battery life by turning off the device after a period of inactivity.
- **Ergonomic Design:** Designed for easy holding and operation.
- **Fever Alert:** Some models provide an audio or visual signal if the measured temperature exceeds a preset fever threshold.

The Echo Dual Function Infrared Thermometer is ideal for home use, pediatric clinics, elder care facilities, schools, or any setting where rapid and hygienic temperature checks are required.

SPECIFICATIONS

- **Display:** HD LED display
- **Measuring range:** Forehead mode: 22°~43°(71.6°F-109.4°F) - Ear mode: 34°~43°(93.2°F-109.4°F) - Object mode: 0°~100°(32°F-212°F)
- **Accuracy:** Forehead mode: $\pm 0.2^{\circ}\text{C}$ (36~39°C); $\pm 0.3^{\circ}\text{C}$ (22~36°C/ 39~43°C); - Ear mode: $\pm 0.2^{\circ}\text{C}$ (36~39°C); $\pm 0.3^{\circ}\text{C}$ (22~36°C/ 39~43°C); - Object mode: $\pm 1.0^{\circ}\text{C}/\pm 2.0^{\circ}\text{F}$
- **Battery:** AAAx2
- **Memory storage:** 20 sets
- **Temperature units:** °C/°F adjustable
- **Automatic shutdown function:** 10s \pm 1s

BROCHURE

[ECHO Thermometer ES-636 Brochure](#)