




ECHO SMART BLOOD PRESSURE MONITOR WITH BLUETOOTH

\$111.00



SKU: ES-625

PRODUCT DESCRIPTION

Smart blood pressure measurement and pulse rate measurement.

FEATURES

- One button operation.
- Measure during inflating technology.
- Irregular heartbeat detection function (IHB).
- Both User 1 and User 2 can store maximum 250 groups of record.
- Guest user option.
- Easy result uploading through Bluetooth.
- TGA Registered (Class II-A) : Australia & New Zealand.
- IOS and Android Compatible.

SPECIFICATIONS

- Power supply: Battery: 6V DC 4x AAA batteries, AC adapter powered mode: 5V 1A (not included)
- Display mode: Digital LCD Display V.A. 60 mm × 64 mm Measurement mode: Oscillographic testing mode Measurement range: Rated cuff pressure: 0 mmHg ~ 299 mmHg Measurement pressure: SYS: 60 mmHg ~ 230 mmHg
- DIA: 40 mmHg ~ 130 mmHg
- Pulse value: (40-199) beat/minute
- Accuracy: Pressure: 5°C ~ 40°C within ± 3 mmHg, pulse value: $\pm 5\%$
- Operating Condition: A temperature range of: +5°C to +40°C

- A relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapour partial pressure greater than 50 hPa
- An atmospheric pressure range of: 700 hPa to 1060 hPa
- Storage & Transportation Condition: Temperature: -20°C to +60°C
- A relative humidity range of $\leq 93\%$, non-condensing, at a water vapour pressure up to 50 hPa
- Measurement perimeter of the upper arm: About 22-32 cm or 22-42 cm