




ELISPOON SELF STABILISING SPOON

\$186.00



SKU: ELI-ELISPOONXL

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

How is ELISpoon different from regular spoons?

ELISpoon is easy to use. Whatever movement of the wrist is performed, the bowl of the spoon remains horizontal and keeps the food in place. The material used is lightweight and strong. A key characteristic is the spoon's affordability. There are no electrical components or batteries, making it carry anywhere and to clean. ELISpoon is a product designed and built to ease the life of the users and their love.

What is the benefit of using ELISpoon?

Using ELISpoon, the individual gains self-confidence and recovers a sense of control. Everyone should enjoy and share moments of conviviality at the dining table, which in essence brings people together and blurs differences. That's what Eli is about. Helping individuals achieve an ability to participate in basic commonplace daily activity.

KEY FEATURES

- ELISpoon is a revolutionary spoon that enables independent eating for people with hand coordination difficulties such as Cerebral Palsy, Spinal Cord Injury, Post-stroke deficits, or Huntington's disease.
- Counterweights and rotation axes keep the spoon bowl level regardless of how your hand or arm twists, bends, or moves, preventing food spillage.
- ELISpoon has no electrical components or batteries, making it easy to use and easy to clean.

Buy this EliSpoon Self Stabilising Spoon online at Medimart as we have NDIS payment method available.

SPECIFICATIONS

Dimensions

Total Length

Total Width

Handle Length

Diameter

Weight

Metal Spoon Head

Colour

ELISpoon XL

Stabilising Soup Spoon

21cm

10.5cm

Approx 12cm

1.9cm

78g

5.5cm (L) x 3.7cm (W)

White