




## ELISPOON SELF STABILISING SPOON



**SKU:** ELI-ELISPOONXL

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

#### How is ELISpoon different from regular spoons?

ELISpoon is easy to use. Whatever movement of the wrist is performed, the bowl of the spoon remains horizontal and keeps the food in place. The material used is lightweight and strong. A key characteristic is the spoon's affordability. There are no electrical components or batteries, making it carry anywhere and to clean. ELISpoon is a product designed and built to ease the life of the users and their love.

#### What is the benefit of using ELISpoon?

Using ELISpoon, the individual gains self-confidence and recovers a sense of control. Everyone should enjoy and share moments of conviviality at the dining table, which in essence brings people together and blurs differences. That's what Eli is about. Helping individuals achieve an ability to participate in basic commonplace daily activity.

### KEY FEATURES

- ELISpoon is a revolutionary spoon that enables independent eating for people with hand coordination difficulties such as Cerebral Palsy, Spinal Cord Injury, Post-stroke deficits, or Huntington's disease.
- Counterweights and rotation axes keep the spoon bowl level regardless of how your hand or arm twists, bends, or moves, preventing food spillage.
- ELISpoon has no electrical components or batteries, making it easy to use and easy to clean.

Buy this EliSpoon Self Stabilising Spoon online at Medimart as we have NDIS payment method available.

## SPECIFICATIONS

<b>Dimensions</b>	<b>ELISpoon XL Stabilising Soup Spoon</b>
<b>Total Length</b>	21cm
<b>Total Width</b>	10.5cm
<b>Handle Length</b>	Approx 12cm
<b>Diameter</b>	1.9cm
<b>Weight</b>	78g
<b>Metal Spoon Head</b>	5.5cm (L) x 3.7cm (W)
<b>Colour</b>	White

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)