



ELK LIFTING CUSHION WITH AIRFLO 24 | EMERGENCY LIFTING CUSHION



SKU: LSS403820

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The ELK Lifting Cushion is a compact, battery powered lifting cushion.

If someone falls over and they can shuffle on to the cushion or be rolled on from the recovery position and assisted to sit upright then the Elk can then be inflated at the touch of a button with someone supporting the individual from behind to bring them to a raise seated position ready to stand.

FEATURES

- Designed to be used with the help of a carer the *ELK* emergency lifting cushion can be used in the tightest of spaces.
- Lifts up to 980 lbs (450 kg)
- Portable weighs just 8lbs and folds up for easy carrying
- Applicable for any age
- · Simple to use
- Minimises moving and handling risk protects carers' backs
- Suitable for indoor and outdoor use
- Airflo controlled (included) The Airflo 24 is a low pressure air compressor which has been designed to inflate the Elk lifting equipment through a 4-Way hand control. The 4-Way hand control is supplied with the Elk.

Click to view the ELK Lifting Cushion with Airflo 24 Brochure.

SPECIFICATIONS

Maximum User Weight: 980lbs Inflated Height: 22"



Inflated Width: 22.5"
Inflated Depth: 20"

Deflated Folded Dimensions: 24.5" x 10" x 4.5"

Weight: 8lbs

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au