

**ENABLE LIFECARE CONFIGURA COMFORT &
ADVANCE ELATA E600C ALTERNATING AIR
CUSHION**



\$2,097.00



SKU: E600C

PRODUCT DESCRIPTION

Purpose & Benefits

The **Elata E600C Alternating Air Cushion** is a dynamic pressure relief solution designed for individuals who have developed pressure injuries. Using alternating inflation and deflation cycles, it redistributes pressure to aid in the prevention and management of skin breakdown.

This cushion is particularly suited to users with pressure injuries in the seating area that have not responded to foam cushions, or those at very high risk due to immobility or the presence of a persistent ulcer. It provides enhanced comfort and accommodates users weighing up to 160 kg.

Features

- Alternating cell inflation and deflation for cyclic pressure redistribution
- Provides comfort and assists with existing pressure injury management
- Compatible with Configura Comfort and Configura Advance chairs

Specifications

- **Includes:** Cushion and pump
- **Maximum User Weight (SWL):** 160 kg

- **Compatibility:** Configura Comfort and Configura Advance chairs